

NexTemp® TraxIt®

Clinical Fever Thermometers

Replaces Digital and Glass Thermometers

NexTemp®

Oral Thermometers



Disposable or Reusable

No Glass

Latex-Free

Mercury-Free

Accurate

**Used in
Hospitals**



TraxIt®

.0 .2 .4 .6 .8
94
95
96
97
98
99
100
101
102
103
104



TraxIt®

Children's Axillary
Thermometers

Wearable

Safe

Birth +

Easy to Use

Reliable

Continuous Reading

Both NexTemp and TraxIt thermometers are made for safe, accurate temperature monitoring and are latex and mercury-free.

MEDICAL
INDICATORS

1589 Reed Rd., Pennington, NJ 08534
888-737-1601 • 609-737-1600 • 609-737-0588 Fax
www.medicalindicators.com

NexTemp® TraxIt®

The fast, accurate, no-hassle thermometers.

Convenient and Economical: NexTemp and TraxIt thermometers are less expensive than any other type of thermometer. They require no special storage conditions, batteries, covers or recalibration processes. NexTemp and TraxIt are available in disposable, reusable or pediatric models. (available in Fahrenheit or Celsius)

Easy To Use: NexTemp and TraxIt thermometers are easier to use and read compared to predecessor dot-matrix thermometers.

Safe: NexTemp and TraxIt replace dangerous glass/mercury, digital electronic and tympanic thermometers. NexTemp and TraxIt thermometers and packaging contain no latex and are non-toxic. They conform to ASTM E1299-96 safety requirements. NexTemp and TraxIt are designed to reduce the risk of cross infection or reinfection that could occur with other thermometer systems.

Accurate: NexTemp and TraxIt color-change thermometers have been proven under strict clinical conditions and meet stringent ASTM accuracy requirements.



NexTemp® Disposable and Reusable Thermometers

Convenient and Versatile: NexTemp can be used for oral, axillary and with the NexTemp sheath, for rectal temperatures.

Directions: Oral Use. Place dot end of the thermometer far back under tongue like a conventional glass thermometer. Wait at least ONE MINUTE before removing. Remove thermometer and READ IMMEDIATELY.

Directions: Axillary Use. Place thermometer high in the axilla (armpit) with the dot-matrix against the torso and parallel to the length of the body. Close arm to hold thermometer in place. Wait at least THREE MINUTES before removing. Remove thermometer and READ IMMEDIATELY.

To Read: The last dot to turn black indicates the correct temperature. The thermometer will reset itself after about 20-30 seconds. Normal oral and axillary temperatures are between 97°F and 99°F (36.1°C and 37.2°C).



TraxIt® Continuous-Reading Wearable Thermometers

Continuous Reading: TraxIt wearable axillary thermometers deliver continuous readings. The first reading is available 2-3 minutes after the thermometer is first applied. Thereafter, readings are virtually instant - just lift the arm to read the temperature as often as you like. Even while the patient is asleep.

Wearable: TraxIt will adhere to the armpit and may be left in place on normal, healthy, intact, hairless skin for up to 48 hours.

Directions for Use: Peel off thermometer from backing and place on the side of the body deep in the axilla, just below the crease in the arm. Close arm and wait 2-3 minutes for the first reading. Lift the arm to read. Thereafter, readings are instant - just lift the arm to read the next temperature. Normal axillary temperature range is 97.0°F to 98.6°F (36.1°C to 37.0°C).

To Read: The last dot to turn BLACK indicates the correct temperature. See example below.

To Remove: Gently peel off and discard.



ORDERING INFORMATION

Disposable Oral/Axillary Thermometers

Individually Wrapped		
NexTemp Single-Use Latex-Free Clinical Thermometers (100 Units/Box)	1112-30°F 2112-20°C	2000 Units (20 Box Case)
Bulk Packaged/Unwrapped		
NexTemp Single-Use Latex-Free Clinical Thermometer (250 Units /Box)	1111-30°F 2111-20°C	5000 Units (20 Box Case)

Reusable Oral/Axillary Thermometers

Packaged in Wallet-Sized Carrier		
NexTemp Latex-Free Personal Thermometer	1211-30°F 2211-20°C	1000 Units (20 Box Case)
Medical Indicators Imprint (50/Box)		
Packed in Sturdy Plastic Case	3072-12°F 3172-12°C	72 Units (6 Box Case)

Rectal Thermometers

NexTemp Plus Single-Use Latex-Free Thermometer with Sheath Cover (100 Units/Box)	1117-30°F 2117-20°C	2000 Units (20 Box Case)
--	------------------------	-----------------------------

Underarm Thermometers

TraxIt Wearable Stick-on Latex-Free Underarm Thermometers (100 Units/Box)	1315-30°F 2315-20°C	2000 Units (20 Box Case)
	1315-14°F 2315-14°C	400 Units (4 Box Case)

Accuracy: All NexTemp and TraxIt thermometers.
±0.2°F (0.1°C) within the critical range of 98.0 to 102.0°F (37.0 to 39.0°C); ±0.4°F (0.2°C) elsewhere.

STORE IN A COOL DRY PLACE



MEDICAL INDICATORS

© 2004 Medical Indicators, Inc.
1589 Reed Rd., Pennington, NJ 08534
Made in U.S.A. All rights reserved.

Questions or comments please contact
Customer Service Toll Free: 888.737.1601
Web: www.medicalindicators.com
Email: customerservice@medicalindicators.com
EU Authorized Rep.: Advena Ltd., PO Box 30, Leominster HR502Q UK

This product is protected under patents #5,676,465,
#6,257,759, #6,284,078 and #6,491,774.

Available in Disposable, Reusable and Pediatric Models