

Tempa•DOT®

**CLINICALLY ACCURATE, SINGLE-USE THERMOMETERS OFFERING
 UNSURPASSED ACCURACY & MAXIMUM INFECTION CONTROL**

PRODUCT SPECIFICATIONS

- Clinical accuracy to ± 0.2 °F/ ± 0.1 °C
- Offers a 2-year shelf life
- No calibrating or cleaning required
- Light-weight and compact
- Should be stored at or below 30°C/86°F
- Available in both Fahrenheit and Celsius
- Does not require batteries, cords, probe covers, or cleaning supplies
- Non-toxic; 100% hygienic; hypoallergenic; glass and mercury- free; not made with natural rubber latex; waterproof; unbreakable
- Temperature reading is permanently locked in

ABOUT TEMPADOT® THERMOMETERS

The Tempa•DOT® line of clinically accurate, single-use thermometers was developed to answer the need for a thermometer that provides highly accurate readings, while also protecting patients and healthcare providers from the spread of germs and infections. All of our Tempa•DOT® thermometers are 100% hygienic and hypoallergenic. And they are manufactured in our FDA and ISO Certified facilities, in a controlled and regulated environment, to ensure their safety.

PRODUCT OPTIONS

We offer Tempa•DOT® thermometers in our standard version, as well as two additional product variations: Tempa•DOT® Sterile, which is the same highly accurate oral/axillary thermometer, that has also been sterilized for instances where infection control is of critical importance; and Tempa•DOT® Plus, which are specifically designed for rectal temperature monitoring.

SAFETY & RELIABILITY

To ensure their safety, we test all of our thermometers utilizing sterility or bioburden testing. The sterility rate on our sterilized thermometers is 0, while the bioburden rate on our standard thermometers is <1 CFU - meaning that our standard thermometers have a population of less than one viable microorganism. It is this level of safety and security our customers have come to rely on for years.

